

DAB Substrate Kit

Catalogue No.: IHC0005

Size: 1200 tests

Storage: Store at 2-8°C for one year except for Reagents A at -20°C, away from light.

Background:

FineTest DAB Substrate (3, 3'-diaminobenzidine) produces a brown reaction product in the presence of peroxidase (HRP) enzyme. Substrate Kit contains all of the reagents necessary toprepare.

REAGENTS:

Reagents A	DAB 20×concentrated solution	3ml	-20°C
Reagents B	H2O2 20×concentrated solution	3ml	2-8°C
Reagents C	TBS 20× concentrated solution	3ml	2-8°C

[•] Avoid storing reagents or working solution in strong direct light.

PREPARATION OF SUBSTRATE WORKING SOLUTION:

• To 1.0 ml of distilled water

Add 50 µl of Reagents C

Add 50 µl of Reagents A

Add 50 µl of Reagents B

• Mix well before use

INSTRUCTIONS FOR USE:

For Tissues or Cells

After incubation with a peroxidase (HRP) detection system, rinse well. Incubate with the substrate working solution at room temperature for 2-10 minutes. Optimal development times should be determined by the investigator.

Wash for 5 minutes in water.

NOTES:

- We recommend using glass-distilled water in the preparation of substrate buffer. Deionized water may con-tain inhibitors of the peroxidase reaction. Solutions containing sodium azide or other inhibitors of peroxidase activity should not be used in diluting the peroxidase substrate.
- DAB are carcinogens. Appropriate care should be exercised when using these reagents including



gloves, eye protection, lab coats, and good laboratory procedures.

For Laboratory Use!