

ACK Lysis Buffer

Catalog No.: K080

Size: 100mL/200mL/500mL

Kit components:

Reagents	100mL	200mL	500mL	Storage
10×ACK Lysis Buffer	100mL	200mL	500mL	2-8°C

Storage:

2-8°C for 12 months

Introduction:

Red Blood Cell Lysis Buffer (RBC Lysis), also known as ACK Lysis Buffer, is a kind of solution used to lyse and remove enucleated red blood cells from blood or tissue samples of human or mouse.

This product is sterile. After treating with this buffer, the single cell suspensions can be used for subsequent primary culture, cell fusion, extraction of nucleic acid or protein, and various routine analysis and detection.

Instructions:

10×ACK Lysis Buffer is concentrated, dilute with ddH₂O to 1×ACK Lysis Working Buffer before use.

Note: 1×ACK Lysis Buffer can store at RT for one week.

Assay Protocol:

For Mouse Spleen or Bone Marrow Samples

1. Centrifuge single cell suspension of the sample at 300 g for 5 min, discard the supernatant.
2. Add 2 mL of 1 ×ACK Lysis Working Buffer to the sample to resuspend the cells, and incubate at RT for 2~3 min (Incubation time should be adjusted accordingly to the actual situation).
3. Add 10 mL PBS to stop the lysis reaction.
4. Centrifuge at 300 g for 5 min, discard the supernatant.
5. Resuspend the cells with PBS or similar physiological buffer for subsequent experiments.

For Human Peripheral Blood Samples

1. Add 100 µL fresh blood into a centrifuge tube, add 2 mL of 1 ×ACK Lysis Working Buffer and mix thoroughly.
2. Put the centrifuge tube at 2-8°C for 10 min, until the fresh blood becomes transparent.

3. Centrifuge at 300 g for 5 min, discard the supernatant. Use a pipette to remove the residual supernatant carefully.
4. Add 2 mL PBS to the tube and repeat step 3.
5. Resuspend the cells with PBS or similar physiological buffer for subsequent experiments.

Note:

1. This product is sterile. Please open and use it in laminar flow bench.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.
3. This reagent is for research use only.