

Intracellular Fixation/Permeabilization Buffer Kit

Catalog No.: K082

Size: 50T/100T/500T

Kit components:

Reagents	50T	100T	500T	Storage
Fixation buffer	10mL	20mL	50mL*2	2-8°C
Permeabilization Buffer (5×)	15mL	30mL	50mL*3	2-8°C

Storage:

2-8°C for 12 months

Introduction:

Intracellular Fixation/Permeabilization Buffer Kit has been formulated and optimized for staining intracellular antigens such as cytokines and chemokines.

Instructions:

Dilute Permeabilization Buffer (5×) with deionized water to 1×Permeabilization Working Solution before use.

Note: It is recommended that 1×permeabilization working solution should be prepared before use and used up within 3 days as far as possible.

Assay Protocol:

1. Take 1×10^6 cells in 100 μ L suspension into the tube per sample.
2. [Optional] Stain cells with a Fixable Viability Dye (self-prepared).
3. [Optional] Block Fc receptors in cell suspensions according to experimental requirements.
4. Stain cell surface markers as need.
5. After incubating with the cell surface marker, add 2 mL of PBS (with 1% BSA, self-prepared) or Cell Staining Buffer [K079], centrifuge at $300 \times g$ for 5 min, discard the supernatant.

6. Resuspend the cells with 200 μ L of PBS (with 1% BSA, self-prepared) or Cell Staining Buffer [K079]. Then add 200 μ L of Fixation Buffer, incubate the cells at room temperature for 30~60 min in the dark (please extend the incubation time to 60 min when the room temperature is lower than 25°C).
7. Add 1 mL of 1 \times Permeabilization Working Solution to each tube and mix fully, centrifuge at 600 \times g for 5 min and discard the supernatant.
8. Resuspend the cells with 100 μ L of 1 \times Permeabilization Working Solution. Add the appropriate volume of intracellular antibody or corresponding isotype control and incubate at least 30 min at room temperature in the dark.
9. Add 2 mL of PBS (with 1% BSA) or Cell Staining Buffer [K079] to each tube and centrifuge at 600 \times g for 5 min, discard the supernatant.
10. Resuspend the cells with appropriate PBS (with 1% BSA) or Cell Staining Buffer [K079], then analyze the samples by flow cytometry.

Note:

- 1.It is normal for the Permeabilization Buffer (5 \times) to have precipitation, and it will not affect the use effect.
- 2.For samples with red blood cells, please lyse red blood cells first.
- 3.The fixation and permeabilization steps may alter the light scatter properties of cells and may increase non-specific background staining. The addition of BSA or fetal calf serum (FBS) in the staining buffer help to reduce non-specific background. It is recommended to use Fixable Viability Dye to eliminate the interference of dead cells in the data analysis process.
- 4.This product is compatible with most commercially available flow antibodies, but some antigenic determinants are sensitive to fixatives and the fixation time needs to be optimized for the situation.
- 5.For your safety and health, please wear the lab coat and disposable gloves before the experiments.