
HRP conjugated Goat Anti-Horse IgG (H+L)

Product Information

Catalog No.:	FNSA-0013
Size:	100µl/500µl
Form:	liquid
Conjugation:	HRP purified from horseradish root
Stabilizer:	3mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.