

anti- ABCB6 antibody

Product Information

Catalog No.:	FNab00031	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

Binds heme and porphyrins and functions in their ATP-dependent uptake into the mitochondria.Plays a crucial role in heme synthesis.

Immunogen information

Immunogen:	ATP-binding cassette, sub-family B(MDR/TAP), member 6
Synonyms:	ABCB6, Mt ABC transporter 3, P glycoprotein related protein, MTABC3, umat
Observed MW:	94 kDa
UniprotID :	Q9NP58


Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB		
Recommended dilution:WB: 1:500-1:5000			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00031(ABCB6 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>