

# anti- ACYP1 antibody

#### **Product Information**

Catalog No.: FNab00133

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Its physiological role is not yet clear.

## **Immunogen information**

Immunogen: acylphosphatase 1, erythrocyte(common) type

Synonyms: ACYP1, ACYPE, Acylphosphatase 1

Observed MW: 11 kDa UniprotID: P07311

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





human brain tissue were subjected to SDS PAGE followed by western blot with FNab00133(ACYP1 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com