

anti- ADAM19 antibody

Product Information

Catalog No.: FNab00140

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Participates in the proteolytic processing of beta-type neuregulin isoforms which are involved in neurogenesis and synaptogenesis, suggesting a regulatory role in glial cell. Also cleaves alpha-2 macroglobulin. May be involved in osteoblast differentiation and/or osteoblast activity in bone(By similarity).

Immunogen information

Immunogen: ADAM metallopeptidase domain 19(meltrin beta)

Synonyms: ADAM 19, ADAM19, FKSG34, MADDAM, Meltrin beta, MLTNB

Observed MW: 105 kDa, observed 70 kDa

UniprotID: Q9H013

Application

Reactivity: Human, Mouse Tested Application: ELISA, WB, IP

Recommended dilution: WB: 1:200-1:1000; IP: 1:200-1:1000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1



Image:



IP Result of anti-ADAM19 (IP:FNab00140, 4ug; Detection:FNab00140 1:300) with NIH/3T3 cells lysate 1200ug.



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab00140(ADAM19 Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com