

anti- ADK antibody

Product Information

Catalog No.:	FNab00172	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

ATP dependent phosphorylation of adenosine and other related nucleoside analogs to monophosphate derivatives. Serves as a potential regulator of concentrations of extracellular adenosine and intracellular adenine nucleotides.

Immunogen information

Immunogen:	adenosine kinase
Synonyms:	ADK, AK, adenosine kinase, Adenosine 5 phosphotransferase
Observed MW:	41 kDa
Uniprot ID :	P55263

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA		
Recommended dilution: None			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com