

anti- AHNAK antibody

Product Information

Catalog No.:	FNab00230
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be required for neuronal cell differentiation.

Immunogen information

Immunogen:	AHNAK nucleoprotein
Synonyms:	AHNAK, AHNAK nucleoprotein, AHNAKRS, Desmoyokin, PM227
Observed MW:	629 kDa
Uniprot ID :	Q09666

Application

Reactivity:	Human		
Tested Application:	ELISA, IF		
Recommended dilution: IF: 1:50-1:500			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889