

# anti- AMOTL1 antibody

## **Product Information**

Catalog No.:	FNab00369	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

### Background

Inhibits the Wnt/beta-catenin signaling pathway, probably by recruiting CTNNB1 to recycling endosomes and hence preventing its translocation to the nucleus.

#### **Immunogen information**

Immunogen:	angiomotin like 1
Synonyms:	AMOTL1, Angiomotin like protein 1, angiomotin like 1, JEAP
Observed MW:	107 kDa
UniprotID :	Q8IY63

# Application


Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB			
Recommended dilution: WB: 1:200-1:2000				
Image:				

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)