

anti- AMPK gamma2 antibody

Product Information

Catalog No.:	FNab00380
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

Immunogen information

Immunogen:	protein kinase, AMP-activated, gamma 2 non-catalytic subunit
Synonyms:	AAKG, AAKG2, AMPK gamma2, AMPK subunit gamma 2, CMH6, H91620p, PRKAG2, WPWS
Observed MW:	63 kDa
UniprotID :	Q9UGJ0

Application


Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, IHC, WB, IF	
Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200		
Image:		

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab00380(PRKAG2 Antibody) at dilution of 1:50

mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00380(PRKAG2 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>