

anti- ANKMY2 antibody

Product Information

Catalog No.: FNab00403

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in the trafficking of signaling proteins to the cilia.

Immunogen information

Immunogen: ankyrin repeat and MYND domain containing 2

Synonyms: ANKMY2, DKFZp564O043, ZMYND20

Observed MW: 49 kDa UniprotID: Q8IV38

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





HeLa cells were subjected to SDS PAGE followed by western blot with FNab00403(ANKMY2 antibody) at dilution of 1:500