

# anti- APOC4 antibody

### **Product Information**

Catalog No.:	FNab00498	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

# Background

May participate in lipoprotein metabolism.

# Immunogen information

Immunogen:	apolipoprotein C-IV
Synonyms:	Apo CIV, ApoC IV, APOC4, apolipoprotein C IV, Apolipoprotein C4
Observed MW:	17 kDa
UniprotID :	P55056

# Application


Reactivity:	Human		
Tested Application:	ELISA, WB, IF		
Recommended dilution: WB: 1:200-1:1000; IF: 1:10-1:100			
Image:			

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com