

# anti- ARMER antibody

### **Product Information**

Catalog No.: FNab00591

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

# Immunogen information

Immunogen: ADP-ribosylation factor-like 6 interacting protein 1

Synonyms: Aip 1, AIP1, ARL 6 interacting protein 1, ARL6IP, ARL6IP1, ARMER,

**KIAA0069** 

Observed MW: 22-24kd UniprotID: Q15041

### **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:1000

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab00591(ARMER antibody) at dilution of 1:300