

# anti- B4GALT7 antibody

### **Product Information**

Catalog No.:	FNab00777	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

### Background

Required for the biosynthesis of the tetrasaccharide linkage region of proteoglycans, especially for small proteoglycans in skin fibroblasts.

## Immunogen information

Immunogen:	xylosylprotein beta 1,4-galactosyltransferase, polypeptide 7(galactosyltransferase I)
Synonyms:	B4GAL T7, B4GALT7, Beta 1,4 GalTase 7, beta4Gal T7, XGALT 1, XGALT1, XGPT, XGPT1
Observed MW:	37kd
Uniprot ID :	Q9UBV7

### Application

Reactivity:	Human	
Tested Application:	ELISA, WB	
Recommended dilution: WB: 1:500-1:1000		

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com



Image:



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00777(B4GALT7 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com