

anti- BAGE4 antibody

Product Information

Catalog No.: FNab00790

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

BAGE4, also named as MLL3P and CT2.4, belongs to the BAGE family. It is a candidate gene encoding tumor antigens. BAGE4 is Cancer/testis antigen 2.4. The antibody can recognize BAGE4, BAGE, BAGE3, BAGE5.

Immunogen information

Immunogen: B melanoma antigen family, member 4

Synonyms: B melanoma antigen 4, BAGE4, Cancer/testis antigen 2.4, CT2.4,

MLL3P

Observed MW: None

Uniprot ID: Q86Y28

Application

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



Immunohistochemical of paraffin-embedded human lung cancer using FNab00790(BAGE4 antibody) at dilution of 1:50

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>