

# anti- BHLHE40 antibody

### **Product Information**

Catalog No.:	FNab00888
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

#### Background

This gene encodes a basic helix-loop-helix protein expressed in various tissues. The encoded protein can interact with ARNTL or compete for E-box binding sites in the promoter of PER1 and repress CLOCK/ARNTL's transactivation of PER1. This gene is believed to be involved in the control of circadian rhythm and cell differentiation.

# **Immunogen** information

Immunogen:	basic helix-loop-helix family, member e40
Synonyms:	BHLHB2, DEC1, SHARP2, STRA13
Observed MW:	48 kDa
Uniprot ID :	O14503

# Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC, IF	
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200		

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889


www.fn-test.com



# Image:



Immunohistochemistry of paraffin-embedded rat brain using FNab00888( BHLHE40 Antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00888( BHLHE40 Antibody) at dilution of 1:1000

# 2

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)