

# anti- BPIL1 antibody

### **Product Information**

Catalog No.:	FNab00937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

# **Immunogen information**

Immunogen:	bactericidal/permeability-increasing protein-like 1
Synonyms:	None
Observed MW:	49kd
Uniprot ID :	Q8N4F0

## Application


Reactivity:	Human			
Tested Application:	ELISA, WB, IHC, IF			
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200				
Image:				

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human lung cancer using FNab00937(BPIL1 antibody) at dilution of 1:100

A549 cells were subjected to SDS PAGE followed by western blot with FNab00937(BPIL1 antibody) at dilution of 1:400



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>