

## anti- Cystatin C antibody

### Product Information

Catalog No.:	FNab02181
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	cystatin C
Synonyms:	ARMD11, CST3, Cystatin 3, cystatin C, CYTC, Gamma trace, Post gamma globulin
Observed MW:	15 kDa
Uniprot ID :	P01034

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab02181(Cystatin C Antibody) at dilution of 1:200



Immunofluorescent analysis of HeLa cells using FNab02181 (Cystatin C antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



A549 cells were subjected to SDS PAGE followed by western blot with FNab02181(CST3 antibody) at dilution of 1:1000