

anti- DUSP22 antibody

Product Information

Catalog No.:	FNab02567
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Activates the Jnk signaling pathway. Dephosphorylates and deactivates p38 and stress-activated protein kinase/c-Jun N-terminal kinase(SAPK/JNK)(By similarity).

Immunogen information

Immunogen:	dual specificity phosphatase 22
Synonyms:	JSP1, LMWDSP2, MKPX
Observed MW:	21 kDa
UniprotID :	Q9NRW4

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human lung cancer using FNab02567(DUSP22 antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02567(DUSP22 Antibody) at dilution of 1:1000