

anti-Fibromodulin antibody

Product Information

Catalog No.:	FNab03121
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7B2
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	fibromodulin
Synonyms:	fibromodulin, FM, FMOD, KSPG fibromodulin, SLRR2E
Observed MW:	43 kDa
UniprotID :	Q06828

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:100 - 1:600
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue using FNab03121(Fibromodulin Antibody) at dilution of 1:300