

## anti- FMNL1 antibody

### Product Information

Catalog No.:	FNab03167
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in the control of cell motility and survival of macrophages(By similarity). Plays a role in the regulation of cell morphology and cytoskeletal organization.Required in the cortical actin filament dynamics and cell shape.

### Immunogen information

Immunogen:	formin-like 1
Synonyms:	C17orf1, C17orf1B, FMNL, FRL1
Observed MW:	150-160 kDa
UniprotID :	O95466

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03167(FMNL1 antibody) at dilution of 1:300