

anti- FUS/TLS antibody

Product Information

Catalog No.:	FNab03242
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	RNA-binding protein FUS
Synonyms:	75 kDa DNA pairing protein, FUS, FUS CHOP, FUS1, FUS-CHOP, hnRNP P2, Oncogene FUS, Oncogene TLS, POMp75, RNA binding protein FUS, TLS, TLS/CHOP
Observed MW:	73 kDa
UniprotID :	P35637

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IP, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab03242(FUS/TLS Antibody) at dilution of 1:200



IP Result of anti-FUS (IP:FNab03242, 3ug; Detection:FNab03242 1:1000) with K-562 cells lysate 4000ug.



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03242(FUS/TLS Antibody) at dilution of 1:2000