

anti- Galectin 9 antibody

Product Information

Catalog No.:	FNab03313
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	lectin, galactoside-binding, soluble, 9
Synonyms:	None
Observed MW:	45kd
UniprotID :	O00182

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-Galectin 9 (IP:FNab03313, 4ug;
Detection:FNab03313 1:500) with SGC-7901 cells
lysate 2400ug.



U-937 cells were subjected to SDS PAGE
followed by western blot with
FNab03313(Galectin 9 antibody at dilution of
1:600