

anti- HNF1B antibody

Product Information

Catalog No.:	FNab03942
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Immunogen information

Immunogen:	HNF1 homeobox B
Synonyms:	FJHN, HNF 1 beta, HNF 1B, HNF1 homeobox B, HNF1B, HNF1beta, HNF2, Homeoprotein LFB3, HPC11, LF B3, LFB3, MODY5, TCF 2, TCF2, Transcription factor 2, VHNF1
Observed MW:	58-60 kDa
Uniprot ID :	P35680

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03942(HNF1B Antibody) at dilution of 1:400



L02 cells were subjected to SDS PAGE followed by western blot with FNab03942(HNF1B antibody) at dilution of 1:800