

anti- MTG1 antibody

Product Information

Catalog No.:	FNab05399
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in the regulation of the mitochondrial ribosome assembly and of translational activity. Displays mitochondrial GTPase activity.

Immunogen information

Immunogen:	mitochondrial GTPase 1 homolog(S. cerevisiae)
Synonyms:	GTPBP7
Observed MW:	37 kDa
UniprotID :	Q9BT17

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05399(MTG1 antibody) at dilution of 1:300