

anti- NPY antibody

Product Information

Catalog No.:	FNab05829
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

Immunogen information

Immunogen:	neuropeptide Y
Synonyms:	CPON, Neuropeptide tyrosine, neuropeptide Y, NPY, PYY4
Observed MW:	11kd
Uniprot ID :	P01303

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05829 (Neuropeptide Y antibody) at dilution of 1:200