

anti- PITX2 antibody

Product Information

Catalog No.:	FNab06474
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Controls cell proliferation in a tissue-specific manner and is involved in morphogenesis. During embryonic development, exerts a role in the expansion of muscle progenitors. May play a role in the proper localization of asymmetric organs such as the heart and stomach. Isoform PTX2C is involved in left-right asymmetry the developing embryo(By similarity).

Immunogen information

Immunogen:	paired-like homeodomain 2
Synonyms:	ARP1, RGS, RIEG, RIEG1
Observed MW:	None
Uniprot ID :	Q99697

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com