

anti-Plakophilin 3 antibody

Product Information

Catalog No.:	FNab06518
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	plakophilin 3
Synonyms:	None
Observed MW:	87 kDa
UniprotID :	Q9Y446

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-Plakophilin 3 (IP:FNab06518, 4ug; Detection:FNab06518 1:600) with A549 cells lysate 1000ug.



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06518(Plakophilin 3 antibody) at dilution of 1:600