

anti- PRPF3 antibody

Product Information

Catalog No.:	FNab06814
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

Immunogen information

Immunogen:	PRP3 pre-mRNA processing factor 3 homolog(<i>S. cerevisiae</i>)
Synonyms:	HPRP3, PRP3
Observed MW:	90 kDa or 78 kDa
UniprotID :	O43395

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-PRPF3 (IP:FNab06814, 5ug;
Detection:FNab06814 1:500) with mouse liver
tissue lysate 4000ug.



Y79 cells were subjected to SDS PAGE followed
by western blot with FNab06814(PRPF3 antibody)
at dilution of 1:1000