

## anti- PTPRO antibody

### Product Information

Catalog No.:	FNab06946
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Possesses tyrosine phosphatase activity. Plays a role in regulating the glomerular pressure/filtration rate relationship through an effect on podocyte structure and function(By similarity).

### Immunogen information

Immunogen:	protein tyrosine phosphatase, receptor type, O
Synonyms:	GLEPP1, PTPU2
Observed MW:	None
Uniprot ID :	Q16827

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)