

---

## anti- RDM1 antibody

### Product Information

Catalog No.:	FNab07221
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May confer resistance to the antitumor agent cisplatin. Binds to DNA and RNA.

### Immunogen information

Immunogen:	RAD52 motif 1
Synonyms:	RAD52B
Observed MW:	32 kDa
UniprotID :	Q8NG50

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis using FNab07221(RDM1 antibody) at dilution of 1:100



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab07221(RDM1 antibody) at dilution of 1:500