

anti- RNF38 antibody

Product Information

Catalog No.:	FNab07356
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as an E3 ubiquitin-protein ligase able to ubiquitinate p53/TP53 which promotes its relocalization to discrete foci associated with PML nuclear bodies. Exhibits preference for UBE2D2 as a E2 enzyme.

Immunogen information

Immunogen:	ring finger protein 38
Synonyms:	None
Observed MW:	43-58 kDa
UniprotID :	Q9H0F5

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07356(RNF38 Antibody) at dilution of 1:300