

## anti- RPB5 antibody

### Product Information

Catalog No.:	FNab07399
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	polymerase(RNA) II(DNA directed) polypeptide E, 25kDa
Synonyms:	None
Observed MW:	25 kDa
Uniprot ID :	P19388

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000;; IHC: 1:50-1:500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer using FNab07399(POLR2E antibody) at dilution of 1:100



IP Result of anti-RPB5 (IP:FNab07399, 4ug; Detection:FNab07399 1:500) with HeLa cells lysate 880ug.

Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07399(POLR2E antibody) at dilution of 1:1000

