
anti- SGT A antibody

Product Information

Catalog No.:	FNab07817
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9E10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Co-chaperone that binds directly to HSC70 and HSP70 and regulates their ATPase activity.

Immunogen information

Immunogen:	small glutamine-rich tetratricopeptiderepeat(TPR)-containing, alpha
Synonyms:	SGT, SGT1
Observed MW:	35 kDa
UniprotID :	O43765

Application


Reactivity:	Human
Tested Application:	ELISA, IF, WB, IHC
Recommended dilution:	WB: 1:1000-1:8000; IF: 1:20-1:200; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab07817(SGTA Antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07817(SGTA Antibody) at dilution of 1:4000