

---

## TOM1 antibody

### Product Information

Catalog No.:	FNab08855
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


May be involved in intracellular trafficking. Probable association with membranes.

### Immunogen information

Immunogen:	target of myb1
Synonyms:	Target of Myb1 membrane trafficking protein, TOM1, Target of Myb protein 1
Observed MW:	54 kDa
Uniprot ID :	O60784

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab08855(TOM1 antibody) at dilution of 1:600