

## anti- VPS72 antibody

### Product Information

Catalog No.:	FNab09450
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Could be a DNA-binding transcriptional regulator. May be involved in chromatin modification and remodeling.

### Immunogen information

Immunogen:	vacuolar protein sorting 72 homolog(S. cerevisiae)
Synonyms:	CFL1, Protein YL 1, Swc2, TCFL1, Transcription factor like 1, VPS72, YL 1, YL1
Observed MW:	42 kDa, observed 57 kDa
UniprotID :	Q15906

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab09450(VPS72 antibody) at dilution of 1:800