

anti- WDR7 antibody

Product Information

Catalog No.:	FNab09501
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	WD repeat domain 7
Synonyms:	KIAA0541, Rabconnectin 3 beta, TRAG, WD repeat containing protein 7, WD repeat domain 7, WDR7
Observed MW:	145 kDa
Uniprot ID :	Q9Y4E6

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis slide using FNab09501(WDR7 Antibody) at dilution of 1:50



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab09501(WDR7 Antibody) at dilution of 1:1000