

## anti- Pan Keratin antibody

### Product Information

Catalog No.:	FNab09935
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells, which are classified into two major sequence types. This antibody is a pan-keratin antibody.

### Immunogen information

Immunogen:	keratin 5
Synonyms:	58 kDa cytokeratin, CK 5, CK5, Cytokeratin 5, K5, keratin 5, KRT5, KRT5A, pan-keratin, Type II keratin Kb5
Observed MW:	46-58 kDa
Uniprot ID :	P13647

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000 - 1:10000; IHC: 1:500 - 1:3000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab09935 (Pan Keratin antibody) at dilution of 1:1000



A431 cells were subjected to SDS PAGE followed by western blot with FNab09935 (Pan Keratin antibody) at dilution of 1:5000