

anti- PMP22 antibody

Product Information

Catalog No.:	FNab10061
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	peripheral myelin protein 22
Synonyms:	PMP22, CIDP, CMT1A, CMT1E, DSS, GAS-3, GAS3, HMSNIA, HNPP, Sp110, peripheral myelin protein 22
Observed MW:	22 kDa
Uniprot ID :	Q01453


Application

Reactivity:	Human,Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500			
Image:			


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Immunohistochemistry of paraffin-embedded human colon cancer using FNab10061(PMP22 antibody) at dilution of 1:100



PC12 cells were subjected to SDS PAGE followed by western blot with FNab10061(PMP22 antibody) at dilution of 1:1000.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com