

anti- RECK antibody

Product Information

Catalog No.:	FNab10074
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Negatively regulates matrix metalloproteinase-9(MMP-9) by suppressing MMP-9 secretion and by direct inhibition of its enzymatic activity. RECK down-regulation by oncogenic signals may facilitate tumor invasion and metastasis. Appears to also regulate MMP-2 and MT1-MMP, which are involved in cancer progression.

Immunogen information

Immunogen:	reversion-inducing-cysteine-rich protein with kazal motifs
Synonyms:	RECK, ST15, reversion-inducing cysteine-rich protein with Kazal motifs
Observed MW:	110 kDa
Uniprot ID :	O95980


Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Image:



Immunohistochemistry of paraffin-embedded human kidney using FNab10074(RECK Antibody) at dilution of 1:100