

anti- CRCM1 antibody

Product Information

Catalog No.:	FNab10170
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	Calcium release-activated calcium channel protein 1
Synonyms:	CRACM1, TMEM142A
Observed MW:	33 kDa
Uniprot ID :	Q96D31

Application

Reactivity:	Rat
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded Rat testis tissue slide using FNab10170(CRCM1 antibody) at dilution of 1:200