

# anti- TSN12 antibody

### **Product Information**

Catalog No.:	FNab10176
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

#### Background

None

## Immunogen information

Immunogen:	Tetraspanin-12
Synonyms:	NET2, TM4SF12
Observed MW:	22 kDa, 30-35 kDa
Uniprot ID :	O95859

## Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200			
Image:			


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10176(TSN12 antibody) at dilution of 1:200



Rat testis tissue were subjected to SDS PAGE followed by western blot with FNab10176(TSN12 antibody) at dilution of 1:1000

## 2

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)