

anti- sumo antibody

Product Information

Catalog No.:	FNab10242
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	Small ubiquitin-related modifier 1
Synonyms:	SMT3C, SMT3H3, UBL1
Observed MW:	11 kDa
Uniprot ID :	P63165

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution	:WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab10242(sumo antibody) at dilution of 1:100



Immunofluorescent analysis of U2OS cells using FNab10242 (sumo antibody) at dilution of 1:50


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com





A549 cells tissue were subjected to SDS PAGE followed by western blot with FNab10242(sumo antibody) at dilution of 1:1000

3

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com