

anti- SIK2 antibody

Product Information

Catalog No.:	FNab10336
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	Serine/threonine-protein kinase SIK2
Synonyms:	SIK2, LOH11CR1I, QIK, SIK-2, SNF1LK2
Observed MW:	104kDa
Uniprot ID :	Q9H0K1

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:500			
Image:			

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275

www.fn-test.com





Immunohistochemistry of paraffin-embedded rat heart tissue slide using FNab10336(SIK2 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab10336(SIK2 antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com