

# anti- FUT2 antibody

#### **Product Information**

Catalog No.: FNab10367

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen: Galactoside alpha-(1,2)-fucosyltransferase 2

Synonyms: FUT2, B12QTL1, SE, SEC2, Se2, sej

Observed MW: 40 kDa Uniprot ID: Q10981

## **Application**

Reactivity: Human, Rat Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


1

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Various lysates were subjected to SDS PAGE followed by western blot with FNab10367(FUT2 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>