

anti- DDX25 antibody

Product Information

Catalog No.:	FNab10413
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

ATP-dependent RNA helicase. Required for mRNA export and translation regulation during spermatid development(By similarity).

Immunogen information

Immunogen:	DEAD(Asp-Glu-Ala-Asp) box polypeptide 25
Synonyms:	DDX25, DEAD box protein 25, GRTH
Observed MW:	50-55 kDa
Uniprot ID :	Q9UHL0

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

293T cells were subjected to SDS PAGE followed by western blot with FNab10413(DDX25 antibody) at dilution of 1:600

