

anti- ADAM28 antibody

Product Information

Catalog No.: FNab10453

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in the adhesive and proteolytic events that occur during lymphocyte emigration or may function in ectodomain shedding of lymphocyte surface target proteins, such as FASL and CD40L. May be involved in sperm maturation.

Immunogen information

Immunogen: ADAM metallopeptidase domain 28

Synonyms: ADAM 28, ADAM28, eMDC II, eMDCII, MDC L, MDC Lm, MDC Ls,

MDCL

Observed MW: 35-90 kDa Uniprot ID: Q9UKQ2

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



A549 cells were subjected to SDS PAGE followed by western blot with FNab10453(ADAM28 antibody) at dilution of 1:800