

G3BP2 antibody

Product Information

Catalog No.:	FNab10613
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


Scaffold protein that plays an essential role in cytoplasmic stress granule formation which acts as a platform for antiviral signaling.

Immunogen information

Immunogen:	GTPase activating protein (SH3 domain) binding protein 2
Synonyms:	G3BP2, KIAA0660
Observed MW:	60 kDa
Uniprot ID :	Q9UN86

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10613(G3BP2 antibody) at dilution of 1:1000