

RBFOX2 antibody

Product Information


Catalog No.:	FNab10679
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	RNA binding protein fox-1 homolog 2
Synonyms:	RNA binding protein fox-1 homolog 2, RFOX2, RNA-binding motif protein 9, RNA-binding protein 9, Repressor of tamoxifen transcriptional activity, FOX2, HRNBP2, RBM9, RTA
Observed MW:	40-60 kDa
Uniprot ID :	O43251

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab10679(RBFOX2 antibody) at dilution of 1:1000